

# **Molegazer**

## **Astronomy meets dermatology**

**Chris Frohmaier - PyData Southampton - 17 Oct 2023**

Dr Chris Frohmaier

- University of Southampton
- Senior Research Fellow in Physics and Astronomy
- Research Software Engineer: Legacy Survey of Space and Time

 cfro.astro@gmail.com

 chrisfrohmaier



# Context

## Astronomy and Dermatology?

Get more data!

Legacy Survey Space Time




Discover exploding stars (supernova)







20 TB/night - I orchestrate with prefect.io




15 PB after 10 year experiment







# Dermatology




CANCER  
RESEARCH  
UK




Science and  
Technology  
Facilities Council




Dr Rubeta Matin PhD FRCP  
Consultant dermatologist  
Senior Lecturer  
Oxford University Hospital Trust





# Context 2


## Medical need



- 60% of melanoma comes from pre-existing naevus (mole)
- >6/day die of melanoma in UK
- Sequential monitoring of patients is labour intensive and not routinely performed
  - Moles are checked by eye and compared to a PDF printout by a consultant!

# The Data


- 50 NHS patients (Oxford cohort)
  - High risk of developing melanoma
- Clinical setting
- Total Body Photography ~3 months collected over 3 years
- 500 patients with single epoch snapshot photography
- >8,000 individual images taken to-date
- All identifying images are removed before sending to Southampton
- Nikon NEF format




**Project Aim:**  
**Predict which naevi will**  
**evolve into melanoma**

**(work-in-progress)**


Aquire Image







Source Detection



Compare History




Predict classification/  
evolution



# The astronomy solution

## Source Extractor



**We could detect moles... but it wasn't a  
scalable solution**



Research  
Software  
Group



Dr Sam Mangam

Senior Research Software Engineer

University of Southampton



Dr Ed Parkinson

Research Software Engineer

University of Southampton



Dr Sam Mangam



Dr Ed Parkinson


**Improve our mole  
detection algorithms  
within the Python  
framework**

**Create a system for  
dermatologists to label  
our data to train future  
models**

# MoleGazer

# Mole Marshal

# The Moleverse




# Molegazer

Detection skin via hue and saturation, then detect background colour.

Use the Laplacian of Gaussian method

- Good for blob detection and edge detection
- Greyscale the image
- Apply a Gaussian blur, then take the second derivative of the pixel intensity map.
- When crosses 0 consider it an edge.





```
from skimage.feature import blob_log
from skimage.color import rgb2gray
from skimage.io import imread
import numpy as np
import matplotlib.pyplot as plt

webb_deep = imread('./Webb_s_first_deep_field.png') #Open JWST deep field
webb = webb_deep[1000:4000,1000:4000] #Small area for demo

webb_gray = rgb2gray(webb) #Convert to Grayscale

#Compute Laplacian of Gaussian
blobs_log = blob_log(webb_gray, max_sigma=30, num_sigma=10, threshold=.1)

#Size is blob is  $\sqrt{2} \times \text{sigma}$  for 2D image
blobs_log[:, 2] = blobs_log[:, 2] * np.sqrt(2)

fig, axes = plt.subplots(1, 1, figsize=(5, 5))
axes.imshow(webb, origin='lower')
for blob in blobs_log:
    y, x, r = blob
    c = plt.Circle((x, y), r, color='palegreen', linewidth=1, fill=False, alpha=0.7)
    axes.add_patch(c)
plt.show()
```



```
from skimage.feature import blob_log
from skimage.color import rgb2gray
from skimage.io import imread
import numpy as np
import matplotlib.pyplot as plt



webb_deep = imread('./Webb_s_first_deep_field.png') #Open JWST deep field
webb = webb_deep[1000:4000,1000:4000] #Small area for demo



webb_gray = rgb2gray(webb) #Convert to Grayscale

#Compute Laplacian of Gaussian
blobs_log = blob_log(webb_gray, max_sigma=30, num_sigma=10, threshold=.1)

#Size is blob is  $\sqrt{2} \cdot \sigma$  for 2D image
blobs_log[:, 2] = blobs_log[:, 2] * np.sqrt(2)

fig, axes = plt.subplots(1, 1, figsize=(5, 5))
axes.imshow(webb, origin='lower')
for blob in blobs_log:
    y, x, r = blob
    c = plt.Circle((x, y), r, color='palegreen', linewidth=1, fill=False, alpha=0.7)
    axes.add_patch(c)
plt.show()
```





r 3.0																			
d 0.22	d 0.25	d 0.28	d 0.28	d 0.29	d 0.29	d 0.3	d 0.31	d 0.31	d 0.32	d 0.33	d 0.34	d 0.34	d 0.34	d 0.34	d 0.35	d 0.35	d 0.35	d 0.36	d 0.36
r 3.0	r 3.0																		
d 0.36	d 0.36	d 0.36	d 0.37	d 0.37	d 0.37	d 0.38	d 0.39	d 0.39	d 0.39	d 0.4	d 0.4								
r 4.67	r 4.67																		
d 0.33	d 0.33	d 0.34	d 0.35	d 0.35	d 0.35	d 0.39	d 0.39	d 0.3	d 0.33	d 0.34	d 0.37	d 0.37	d 0.39	d 0.3	d 0.32				
r 16.33																			
d 0.35	d 0.39	d 0.31	d 0.32	d 0.34	d 0.37														

surroundings\_difference\_min=0.2

r 3.0	r 3.0	r 3.0	r 3.0	r 3.0	r 3.0	r 3.0	r 3.0	r 3.0	r 3.0	r 3.0	r 3.0	r 3.0	r 3.0	r 3.0	r 3.0	r 3.0	r 3.0	r 3.0	
d 0.4	d 0.4	d 0.4	d 0.4	d 0.41	d 0.42	d 0.42	d 0.43	d 0.43	d 0.43	d 0.43									
r 3.0	r 3.0	r 3.0	r 3.0	r 3.0	r 3.0	r 3.0	r 3.0	r 3.0	r 3.0	r 3.0	r 3.0	r 3.0	r 3.0	r 3.0	r 3.0	r 3.0	r 3.0	r 3.0	r 3.0
d 0.43	d 0.44	d 0.44	d 0.44	d 0.44	d 0.44	d 0.45	d 0.45	d 0.45	d 0.45	d 0.46	d 0.46	d 0.46	d 0.46	d 0.47	d 0.47	d 0.48	d 0.48	d 0.48	d 0.48
r 3.0	r 3.0	r 3.0	r 3.0	r 3.0	r 3.0	r 3.0	r 3.0	r 3.0	r 3.0	r 3.0	r 3.0	r 3.0	r 3.0	r 3.0	r 3.0	r 3.0	r 3.0	r 3.0	r 3.0
d 0.48	d 0.49	d 0.49	d 0.49	d 0.49	d 0.5	d 0.51	d 0.51	d 0.51	d 0.51	d 0.52	d 0.52	d 0.53	d 0.53	d 0.53	d 0.53				
r 3.0	r 3.0	r 3.0	r 3.0	r 3.0	r 3.0	r 3.0	r 3.0	r 3.0	r 3.0	r 3.0	r 3.0	r 3.0	r 3.0	r 3.0	r 3.0	r 3.0	r 3.0	r 3.0	r 3.0
d 0.53	d 0.54	d 0.54	d 0.56	d 0.57	d 0.57	d 0.59	d 0.4	d 0.43	d 0.44	d 0.45	d 0.45	d 0.46	d 0.46	d 0.57	d 0.44	d 0.52	d 0.41	d 0.43	d 0.51
r 16.33	r 16.33																		
d 0.51	d 0.52																		


surroundings\_difference\_min=0.4

r 3.0	r 3.0	r 3.0	r 4.67																
d 0.61	d 0.62	d 0.62	d 0.63	d 0.65	d 0.66	d 0.67	d 0.78	d 0.63	d 0.64	d 0.66	d 0.68	d 0.7	d 0.71	d 0.71	d 0.73	d 0.79	d 0.8	d 0.61	
r 6.33	r 11.33																		
d 0.63	d 0.66	d 0.72	d 0.74	d 0.77	d 0.79	d 0.8	d 0.6												
r 3.0	r 4.67																		
d 0.81	d 0.81	d 0.84	d 0.85	d 0.89	d 0.9	d 0.9	d 0.92	d 0.81	d 0.82	d 0.85	d 0.9	d 0.91	d 0.84	d 0.88					

surroundings\_difference\_min=0.5

r 3.0	r 3.0	r 3.0	r 3.0	r 3.0	r 3.0	r 3.0	r 3.0	r 3.0	r 3.0	r 3.0	r 3.0	r 3.0	r 3.0	r 3.0					
d 0.81	d 0.81	d 0.84	d 0.85	d 0.89	d 0.9	d 0.9	d 0.92	d 0.81	d 0.82	d 0.85	d 0.9	d 0.91	d 0.84	d 0.88					

surroundings\_difference\_min=0



# MoleMarshal

## Labelling the data



- Citizen Science platform engaging the public with research
- Made famous for asking public to label images of galaxies
- MoleMarshal is a private project
  - Patient sensitive data
  - Only accessed from behind UoS firewall
- Consultant dermatologists can look through the data and label it
- (I am not allowed to demo this - patient confidentiality)

# Interface

The screenshot displays the MoleMarshal-Development web application interface. At the top left is a user icon with a circular arrow symbol. Next to it is the text "MoleMarshal-Development". Along the top right are navigation links: ABOUT, CLASSIFY (which is underlined), TALK, COLLECT, RECENTS, and LAB.

The main content area features a large image of a skin lesion with a red circle highlighting a specific area. Below this image is a smaller row of three smaller images, each with a red circle, likely showing different views or stages of the mole.

To the right of the images is a task panel. It includes tabs for "TASK" and "TUTORIAL", with "TASK" currently selected. The task question "Is there a mole in the top image?" has two options: "No" and "Yes". A blue rectangular box labeled "NEED SOME HELP WITH THIS TASK?" contains these buttons. At the bottom of the panel are three buttons: "Done & Talk" (blue), "Done" (green), and a settings gear icon.


At the very bottom center of the page are small icons for information, heart, and list.

# Interface



MoleMarshal-Development

ABOUT CLASSIFY TALK COLLECT RECENTS LAB



Is there a mole in the top image?

**TASK** **TUTORIAL**

No

Yes

NEED SOME HELP WITH THIS TASK?

Done & Talk Done ⚙


SWITCH TO DARK THEME

# Interface



MoleMarshal-Development

ABOUT CLASSIFY TALK COLLECT RECENTS LAB



How severe is this mole, in comparison to the moles in the lower images?

**TASK**

**TUTORIAL**

1

2

3

4


5

N/A


NEED SOME HELP WITH THIS TASK?


Back Done & Talk Done ⚙


SWITCH TO DARK THEME



**Finding similar moles...  
efficiently label data**







Only use the values at these nodes  
Similar regions in this nD-space would  
be occupied by similar looking moles?



```
# for loading/processing the images
from keras.preprocessing.image import load_img
from keras.preprocessing.image import img_to_array
from keras.applications.inception_resnet_v2 import preprocess_input

# models
from keras.applications.inception_resnet_v2 import InceptionResNetV2
from keras.models import Model

# clustering and dimension reduction
from sklearn.cluster import KMeans
from sklearn.decomposition import PCA

# for everything else
import os
import numpy as np
import matplotlib.pyplot as plt
from random import randint
import pandas as pd
import pickle
import glob
```



```
# load model
model = InceptionResNetV2(
    weights='imagenet', # Load weights pre-trained on ImageNet.
    input_shape=(299, 299, 3),
)
model = Model(inputs=model.inputs, outputs=model.layers[-2].output)
```

```
imageList = glob.glob('/my/very/private/medical/data/stamps/*.png')

def extract_features(file, model):

    img = load_img(file, target_size=(299,299))
    img = np.array(img)
    reshaped_img = img.reshape(1,299,299,3)
    imgx = preprocess_input(reshaped_img)

    features = model.predict(imgx, use_multiprocessing=True)

    return features

data = {}
for mole in imageList:
    feat = extract_features(mole,model)
    data[mole] = feat
```




```
pca = PCA(n_components=100, random_state=42)
pca.fit(feat)
x = pca.transform(feat)

numberOfClasses = 10 #Guess how many classes

kmeans = KMeans(n_clusters=numberOfClasses, random_state=42)
kmeans.fit(x)
```

# Plot the clusters

Randomly sampling from each cluster



A more efficient labelling system

# (For fun)

## What did the classifier predict?


tick, 0.91



tick, 0.96



tick, 0.28



banana, 0.18



mashed\_potato, 0.05



# Future work

- Things seem to work on single snapshots!
  - We can detect moles on patients
- More diverse cohort!
- Need to build in time series
  - Tracking moles over multiple years
  - How to align patients so we know we are looking at the same moles
  - Astronomy techniques don't work reliably for this problem.

# Questions?



Dr Chris Frohmaier

 chrisfrohmaier



cfro.astro@gmail.com

 cfrohmaier

